

PTFE Peelable Introducer



Greatbatch Medical offers a broad portfolio of peelable introducers in a wide variety of sizes. Our design features facilitate access in a variety of clinical applications.

DEPENDABLE DESIGN

- Low crack-force and reliable, gliding peel performance aid in maintaining the position of the implanted device
- Smooth sheath-to-dilator transition, resulting in low insertion force and atraumatic vessel access
- Features a low-profile handle and a locking hub system, which prevents dilator pop-out during insertion

HIGHLY EFFICIENT

- French size indicated on dilator hub
- Configurations available in bulk (non-sterile) or procedure kits (sterile)

Greatbatch Medical
10000 Wehrle Drive
Clarence, NY 14031
tel 716-759-5600
fax 716-759-2562

sales@greatbatchmedical.com
www.greatbatchmedical.com

1002672-000 Rev. A
©2010 Greatbatch Medical. All rights reserved.


Intended for use in the percutaneous insertion of pacing leads or catheters in the venous system.

PTFE Peelable Introducer (.018in guidewire compatible)										
Introducer French Size	3.5	3.5	4.5	4.5	5.0	5.0	6.0	6.0	7.0	7.0
Sheath I.D.	.050in	.050in	.066in	.066in	.073in	.073in	.084in	.084in	.097in	.097in
Sheath Length	6cm	11cm	6cm	11cm	6cm	11cm	6cm	11cm	6cm	11cm
Bulk Model Number	10613-001	10613-002	10613-003	10613-004	10613-005	10613-006	10613-007	10613-008	10613-009	10613-010

PTFE Peelable Introducer (.035in/.038in guidewire compatible)									
Introducer French Size	7.0	8.0	9.0	10.0	12.0	13.0	14.0	15.0	16.0
Sheath I.D.	.099in	.112in	.119in	.138in	.159in	.172in	.190in	.198in	.209in
Sheath Length	15cm	15cm	15cm	15cm	15cm	15cm	14cm	14cm	14cm
Bulk Model Number	10491-004	10491-005	10491-006	10491-007	10491-008	10491-009	10491-010	10491-011	10491-012

PTFE Kit Includes:

- Introducer
- 0.035 access wire
- 18 gauge needle
- Syringe



Greatbatch Medical
 10000 Wehrle Drive
 Clarence, NY 14031
 tel 716-759-5600
 fax 716-759-2562

sales@greatbatchmedical.com
www.greatbatchmedical.com