

ALCath Flutter Flux G eXtra

3.5 mm GoldTip ablation catheter



Product Highlights

- **Refined Control**
Superior high-density braiding for direct 1:1 motion transfer
- **Reinforced Stability**
Lateral stabilization of the tip segment for effective tip-to-tissue contact
- **Reduced Length**
Catheter length optimized for the right atrium for enhanced high-precision steering

Ordering Information

Model	Tip	Curve	Spacing	French	Length	Order number
ALCath Flutter Flux Black G eXtra	3.5 mm	Black	2-5-2 mm	7F	95 cm	380783


AlCath Flutter Flux G eXtra

Technical Data

Catheter specification	
Tip electrode size	3.5 mm
Tip electrode type	Irrigated
Tip temperature sensor	Thermocouple
Ring electrode size	1.55 mm
Number of electrodes	4
Electrode spacing	2-5-2 mm
Electrode material	Gold (Au)
Shaft diameter	7F (2.3 mm)
Braiding	High density
Curve type	Black
Curve shape	Steerable – FullCircle
Reach at 90° deflection	73 mm
Length of the deflectable area	95 mm
Effective catheter length	95 cm
Connection plug	9-pole, REDEL Luer-lock (female)
Patient cable	PK-109, 111, 112, 142, 143, 185

Patient cable	Ablation system	Order number
PK-109	Stockert EP-Shuttle	330 079
PK-111	Osypka HAT 300 Smart	330 080
PK-112	Medtronic Atakr, Radionics RFG 3-E	330 081
PK-142	IBI 1500	362 442
PK-143	Boston Scientific Maestro 3000	371 597
PK-185	BIOTRONIK AbControl MDS	345 954

Curve chart	Length L	Reach R
Black	95 mm	73 mm



ALCath Flutter LT G

8mm GoldTip ablation catheter



Product Highlights

- **Refined Control**
Superior high-density braiding for direct 1:1 motion transfer
- **Reinforced Stability**
Lateral stabilization of the tip segment for effective tip-to-tissue contact
- **Reduced Length**
Catheter length optimized for the right atrium for enhanced high-precision steering

Ordering Information

Model	Tip	Curve	Spacing	French	Length	Order number
ALCath Flutter Black LT G	8 mm	Black	2-5-2 mm	7F	95 cm	368997

AlCath Flutter LT G

Technical Data

Catheter specification	
Tip electrode size	8 mm
Tip electrode type	Non-irrigated
Tip temperature sensor	Thermocouple
Ring electrode size	1.55 mm
Number of electrodes	4
Electrode spacing	2-5-2 mm
Electrode material	Gold (Au)
Shaft diameter	7F (2.3 mm)
Braiding	High density
Curve type	Black
Curve shape	Steerable – FullCircle
Reach at 90° deflection	73 mm
Length of the deflectable area	95 mm
Effective catheter length	95 cm
Connection plug	9-pole, REDEL
Patient cable	PK-109, 111, 112, 142, 143, 185

Patient cable	Ablation system	Order number
PK-109	Stockert EP-Shuttle	330 079
PK-111	Osyka HAT 300 Smart	330 080
PK-112	Medtronic Atakr, Radionics RFG 3-E	330 081
PK-142	IBI 1500	362 442
PK-143	Boston Scientific Maestro 3000	371 597
PK-185	BIOTRONIK AbControl MDS	345 954

Curve chart	Length L	Reach R
Black	95 mm	73 mm

